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Cat: G2950

Size: 100mL

Storage: RT, avoid light, valid for 6 months.

Introduction

Turk's Solution, also known as Tuerk's solution, which is used for leukocyte count. The principle is that Turk's solution as a blood diluent, Mature red blood cells are all dissolved. Fill them into the counting chamber, count the number of leukocytes in a certain volume under the microscope, and calculate the number of leukocytes in each liter by conversion.

Turk's Solution is only used in scientific research, not in clinical diagnosis.

Self Provided Materials

Fresh whole blood, Micropipette, Cell counting plate, Microscope

Protocol(for reference only)

1. Take a small EP tube and add 0.38mL Turk's Solution.
2. Take 20 fresh or anticoagulant blood samples with a micro pipette μ l. Wipe off the remaining blood outside the tube, add it to the bottom of Turk solution, discharge the blood, gently suck the supernatant, clean the straw for 2 ~ 3 times, and mix it immediately. The dilution ratio is 1:20. If there are few nucleated cells, it can also be diluted according to 1:10.
3. Let it stand at room temperature for 3 ~ 5min. After mixing again, take 10ul and flush it into the counting cell. Pay attention to prevent bubbles or overflow. Let it stand for 2-5min and count after leukocyte precipitation.
4. Count the number of leukocytes in four corners and four squares in turn under low power microscope. Pressure line cells are counted according to the principle of "counting up and not counting down, counting left and not counting right".

Calculate: leukocyte number /L=(N÷4)×10×20×10⁶= N × 50×10⁶

N	Total number of leukocytes in 4 squares
÷4	Average number of leukocytes in each large square (0.1 μ L)
×10	The volume of a square lattice is 0.1 μ L, which is converted into 1.0 μ L
×20	Blood dilution times
×10 ⁶	Conversion from 1 μ L to 1 L

Note

1. When taking blood, avoid squeezing too much and keep in appropriate cupuncture depth.
2. The small test tube and counting plate shall be clean to avoid misidentification of cells.
3. Within the reference range, the number of cells between squares should not differ by more than 8, and the difference between two repeated counts should not exceed 10%.
4. For your safety and health, please wear experimental clothes and disposable gloves.