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**Cat:** G3272

**Size:** 2×50mL/2×100mL/2×500mL

**Storage:** RT, avoid light, valid for 1 year.

### Kit Components

Reagent	2×50mL	2×100mL	2×500mL	Storage
Reagent A: Weil Iron Hematoxylin Solution A	50mL	100mL	500mL	RT, avoid light
Reagent B: Weil Iron Hematoxylin Solution B	50mL	100mL	500mL	RT, avoid light
Mix A and B in equal to form Weil Hematoxylin Working Solution. It is not suitable to prepare in advance.				

### Introduction

Myelin sheath is a layer of membrane wrapped around the axons of nerve cells, that is, myelin sheath is composed of myelin cells and cell membrane. It is a multilayer lipid double-layer structure formed by the plasma membrane of nerve membrane cells spirally winding along the axis of axon. There is a Langfei's knot on the myelin sheath, which can make the spirit jump and transmit through impulse.

Weil's Iron Hematoxylin Stain Kit can show whether the myelin sheath is complete, denatured, necrotic and repaired under pathological conditions. It has significance for pathological diagnosis and research of nerve tissue. For example, when the nerve fiber is damaged, the myelin sheath may swell, zigzag into a sphere, break or completely disappear without sheath.

*Note: Alum solution, differentiation solution and bluing solution should be prepared by ourselves. For the first operation, we recommend G3270 kit, which has been equipped with all dyeing reagents, saving your time.*

### Protocols(for reference only)

1. Cut into paraffin section, dewax to distilled water.
2. Soak the section in Weil Iron Hematoxylin Solution in a warm box or water bath at 56 °C for 20 min. Rinse with distilled water.
3. Add Iron Alum Solution dropwise for differentiation, while controlling the differentiation time under the microscope (distinguishing between normal myelin sheath and gray matter or degenerative areas), and rinse with distilled water.
4. Add Weil Differentiation Solution (eliminate the background color). Rinse with distilled water.
5. Add Weil Bluing Solution and fully wash with water.
6. Conventionally dehydrate, transparent by xylene and seal with resinene.

### Result

Myelin Sheath	Black
Background	Colorless

### Note

1. This kit is simple and fast. Iron Alum Solution differentiation is a key step. It is necessary to observe the degree of differentiation under the microscope.
2. The best fixative is 10% formalin.
3. When dyeing in a warm box or water bath, pay attention to prevent the evaporation of the dye solution.
4. For your safety and health, please wear experimental clothes and disposable gloves.