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Cat: G3320

Size: 2×50mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent		2×50mL	Storage
Reagent(A):Lillie Stain Solution	A1:Lillie's Ferrous Iron Solution A	25mL	RT
	A2:Lillie's Ferrous Iron Solution B	25mL	RT, avoid light
Mix A1,A2 in equal to form Lillie Stain Solution before use, which can not keep for long.			
Reagent(B): Nuclear Fast Red Solution		50mL	RT, avoid light

Introduction

The Prussian Blue Reaction is a method of staining non heme iron (containing heme, iron transporters, iron deposition, etc.) in cells or tissues containing granular iron. It plays an important role in the validation of diseases such as ischemic anemia, hemolytic anemia, hemochromatosis, and abnormal hemoglobin metabolism, and can also be used in conjunction with other enhancers for the detection of iron death models.

The Prussian blue reaction was first proposed by Perls and is therefore most commonly referred to as the Perl iron staining method. However, since its proposal in 1867, numerous chemists and pathologists have optimized this reaction, among which the Lillie method, which is more widely used and well optimized, can label divalent and trivalent iron separately. This kit is a ready to use Lillie method ferrous iron staining kit, prepared with high-quality raw materials, with stable reagents and good staining specificity.

Protocol(for reference only)

1. The organization recommends fixation in neutral formalin fixative (10%) or other non acidic fixative, conventional dehydration and embedding. (See Note 1)
2. Cut the paraffin section in 5µm thickness. Conventional dewaxing to water. The recommended thickness for frozen sections is 10µm. Regular rewarming and rehydration. (See Note 2)
3. Stain with Lillie Stain Solution for 25-30min and rinsed with distilled water for 2-5min.
4. Re-dye with Nuclear Fast Red for 5-10min. Rinse with distilled water for 1-5s.(See Note 3)
5. Conventionally dehydrate and transparent, then seal with resinene.

Result

Ferric Iron	Blue or Bluish Green
Background	Red

Negative Control(Optional)

Take the same control section and dewax it to water, then treat it with a 1% hydrochloric acid solution at room temperature for 20min before staining normally. The result is negative for staining.

Take some ferrous salt and prepare an aqueous solution samples (commonly 0.1% ferrous sulfate aqueous solution), and drop a drop of staining solution, the result is positive for kit.

Take some iron salt and prepare an aqueous solution samples (commonly 0.1% ferric chloride or ferric sulfate aqueous solution). Drop a drop of the staining solution, the result is a negative for kit.

Note

1. This staining method for tissue sampling should avoid using acidic fixative or chelating agent treatment that may cause iron ion loss. It is recommended to use G2161-Neutral Buffered Formalin Fixative, 10% or P1110-Paraformaldehyde,4% for tissue fixation.
2. This staining method is applicable to paraffin sections, frozen sections, and resin sections. However, due to the instability of divalent iron, it is necessary to avoid contact with oxidants or prolonged treatment during the production process. Therefore, it is usually recommended to prepare sections after tissue perfusion staining. If perfusion is not possible, it is recommended to prepare frozen sections for staining.
3. Reagent (B): Nuclear Fast Red Solution is a colloidal solution.Long term storage at low temperatures (below 25 °C) will produce a small amount of precipitation, which is a normal phenomenon and generally does not affect its use. If obvious turbidity is observed, it can be restored to use by boiling water for 5-10 min.
4. During the entire operation process, it is necessary to avoid iron ion pollution, and distilled water is preferred for cleaning, as tap water often contains iron.
5. For your safety and health, please wear laboratory clothes and disposable gloves.