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**Cat:** G3500

**Size:** 4×50mL/4×100mL

**Storage:** RT, avoid light, valid for 6 months.

### Kit Components

Reagent		4×50mL	4×100mL	Storage
Reagent(A) : Weigert Iron Hematoxylin Solution	A1: Weigert Solution A	25mL	50mL	RT, avoid light
	A2: Weigert Solution B	25mL	50mL	RT, avoid light
Before use, mix A1 with A2 in equal amount to form Weigert Iron Hematoxylin Solution. It is not suitable to prepare in advance.				
Reagent(B): Acidic Differentiation Solution		50mL	100mL	RT
Reagent(C): Pollak Staining Solution		50mL	100mL	RT, avoid light
Reagent(D): Pollak Differentiation Solution		50mL	100mL	RT

### Introduction

Masson Trichrome Staining, also known as Masson staining, is the most classical method in connective tissue staining. It is an authoritative and classic technical method of collagen fiber staining. Masson Trichrome Staining usually refers to the staining of nucleus and selective display of collagen and muscle fibers. The dyeing principle of this method is related to the size of anionic dye molecules and the permeability of tissues: the size of molecules is reflected by molecular weight, and small molecular weight is easy to penetrate the tissues with dense structure and low permeability; while large molecular weight can only enter the tissues with loose structure and high permeability. Pollak Trichrome Staining is a polychromatic staining method of connective tissue improved on Masson Trichrome Staining method. It can stain various components of connective tissue by using mordant and promoter at the same time.

### Protocol(for reference only)

1. For paraffin section, conventionally dewax to distilled water.
2. Stain with prepared Weigert Iron Hematoxylin Solution for 15-20 mins.
3. Differentiate in Acidic Differentiation Solution for 5-10 s.
4. Blue in tap water for 10mins, and wash with distilled water for 2-4 times.
5. Drop with Pollak Staining Solution for 3-7mins. Quickly wash with distilled water.
6. Differentiate by 0.2% glacial acetic acid(control under the microscope) about 10-20s. Quickly wash with distilled water.
7. Quickly dehydrate in ethanol and each level for 3-5s. Transparent by xylene for 2 times, 1-2min each.
8. Seal with resinene.

### Result

Collagen Fiber, Mucus, Cartilage, Nerve Fiber	Blue
Muscle, Elastic Fiber	Red
Cellulose	Purplish Red
Red Blood Cell	Orange Red
Nucleus	Bluish Black

### Note

1. Section dewaxing should be as clean as possible. Weigert Iron Hematoxylin Solution is ready to use, and generally loses its dyeing power in 6 h.
2. Tissue fixation plays an important role. For different fixatives, can prolong or shorten the staining time.
3. The dyeing time of Pollak Staining Solution should be strictly controlled: too short dyeing time leads to too deep red; too long dyeing time leads to too deep green or blue. Over differentiation will make the color of mucus, collagen fiber and cellulose too light to observe.
4. Pollak Differentiation Solution can make the color more clear and bright. If the use is large, G2940- Acetic Acid Solution can be prepared to replace it.
5. For your safety and health, please wear experimental clothes and disposable gloves.