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Cat: G3525

Size: 6×50mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	6×50mL	Storage
Reagent(A):Gordon-Sweets Oxiadant	A1: GS Oxiadant A	25mL RT, avoid light
	A2: GS Oxiadant B	25mL RT
Mix equal parts of GS Oxidant A and B to form Gordon-Sweets Oxiadant before use. It is ready to use.		
Reagent(B): Oxalic Acid Solution	50mL	RT
Reagent(C): Ferric Ammonium Sulfate Solution	50mL	RT, avoid light
Reagent(D):Gordon-Sweets Ammoniacal Silver Solution	50mL	2-8°C, avoid light
Reagent(E): Gordon-Sweets Reductant	50mL	RT
Reagent(F): Nuclear Fast Red Solution	50mL	RT, avoid light

Introduction

Reticular fiber is a kind of fiber in reticular connective tissue, which is produced by reticular cells. Its diameter is mostly 0.2-1.0 μ m. it has toughness but no elasticity. There are many dyeing methods for reticular fiber, but the dyeing principle is basically the same, and most of them adopt ammonia silver staining method. In the traditional method, after reduction, gold chloride solution is used for color matching, and then hypo solution is used to wash away the unreduced silver salt on the tissue. This improved method omits this step and makes the contrast of reticular fiber clearer.

Reticulin Stain Kit(Gordon-Sweet Sliver Method) uses acid oxidant and nuclear fast red re-dyeing solution to get better result. Frozen section, low temperature section and collodion section can be used for reticular fiber staining.

Protocol(for reference only)

Reagent(F): Nuclear Fast Red Solution may produce suspended solids or a small amount of precipitation due to flocculation. It is recommended to take supernatant or boil water bath for 5-10min and then air it to 30-40 °C. (see Note 2)

1. Fix the tissue in 10% formalin fixative, then conventionally dehydrate and embed.
2. Cut into paraffin section in 4 μ m thick, and dewax to distilled water.
3. Drop in the prepared Gordon-Sweets Oxiadant and oxidize for 3min. Rinse with running water for 20s.
4. Add Oxalic Acid Solution for 1min, wash with distilled water twice for 1min each.
5. Add Ferric Ammonium Sulfate Solution for 10min, wash with distilled water twice for 1min each.
6. Add Gordon-Sweets Ammoniacal Silver Solution for 11s, wash with distilled water twice for 2min each.
7. Add Gordon-Sweets Reductant for 2min and rinse in running water for 1min.
8. Add Nuclear Fast Red Solution for 10min and wash slightly with water for 20s.
9. Conventionally dehydrate and transparent, seal with resinene.

Result

Reticular Fiber	Black
Collagen Fiber	Yellow to Yellowish Brown
Nucleus	Red
Cytoplasm	Light Yellow

Note

1. Glass container must be soaked in washing solution for one day, washed with tap water and washed twice with distilled water.
2. Reagent(F): Nuclear Fast Red Solution is a colloidal solution, which is stored at low temperature (lower than 25 °C) or stored for a long time. Suspended solids or a small amount of precipitation are generated due to flocculation, which is a normal phenomenon and generally does not affect the use. If the colloid solution is evenly dispersed in the boiling bath, tighten the bottle cap for 5-10min to recover the turbid solution.
3. 10% formalin fixative is a suitable fixative. It is not suitable to use fixative containing mercury such as Zenker's Solution, otherwise it is easy to cause nonspecific deposition of sections.
4. Gordon-Sweets Ammoniacal Silver Solution should be kept away from light at 4 °C and used in RT.
5. For your safety and health, please wear experimental clothes and disposable gloves.