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Cat: G3550

Size: 6×50mL/6×100mL

Storage: RT, avoid light, valid for 1 year.

Kit Components

Reagent		6×50mL	6×100mL	Storage
Reagent(A): Weigert Iron Hematoxylin Solution	A1: Weigert Solution A	25mL	50mL	RT, avoid light
	A2: Weigert Solution B	25mL	50mL	RT
Mix equal parts of A1 and A2 to form Weigert Iron Hematoxylin Solution before use.				
Reagent(B): Acid Differentiation Solution		50mL	100mL	RT
Reagent(C): Acid Ponceau Solution		50mL	100mL	RT, avoid light
Reagent(D): Weak Acid Solution		50mL	100mL	RT
Reagent(E): Orange G Solution		50mL	100mL	RT, avoid light
Reagent(F): Light Green Solution		50mL	100mL	RT, avoid light

Introduction

There are many methods of bone staining, such as Toluidine Blue method, Alician Blue method, Saffron O method, Goldner Trichrome method and so on. Goldner Trichrome Staining is similar to Masson Trichrome Staining. Trichrome staining usually refers to staining the nucleus and selectively displaying collagen and muscle fibers. This kit is mainly used for staining bone substances and has good effect on cell staining, especially for the study of metabolic diseases (Paget disease, osteodystrophy, etc.), so as to evaluate the activity of osteoblasts and osteoclasts, and easily identify the metastasis tumor cells in bone marrow.

Self Prepared Materials

Series of ethanol, Distilled water

Protocols (for reference only)

1. Slices are dewaxed with xylene or plasticized with organic solvents, and rehydrated with gradient ethanol.
2. Dye with prepared Weigert Iron Hematoxylin Solution for 15-20 min and rinse with running water for 1 min.
3. Quickly differentiate with Reagent(B): Acid Differentiation Solution (generally <5 s). Wash with distilled water for 1 min.
4. Add Acid Ponceau Solution for 5 min.
5. During the above process, prepare a weak acid working solution in a ratio of 4:1 distilled water to Weak Acid Solution, and wash with the weak acid working solution for 15- 30 s.
6. Add Orange G Solution til the collagen fibers become red, usually taking 3-10 min. (See Note 5)
7. Rinse with prepared weak acid working solution for 15- 30 s.
8. Dye directly in Reagent(F): Light Green Solution for 5 min, and rinse 3 times with prepared weak acid working solution for 15 s each time.
9. Absorb or air dry, dehydrate with anhydrous ethanol, and seal with resinene.

Result

Neonatal bone	Orange-Red
Osteoid	Purple
Mature bone	Green
Nucleus	Blue-Grey

Note

1. Section dewaxing should be as clean as possible. Fixation plays an important role. According to different fixatives, can prolong or shorten the dyeing time.
2. The Weigert Iron Hematoxylin Solution is ready to use, which generally loses its dyeing power for 24 h.
3. The differentiation time of Acid Differentiation Solution should be determined according to the thickness of slice, the type of tissue and the old and new.
4. Weak Acid Solution can make the color more clear and bright. If the use amount is large, can replace with 0.1-1% acetic acid solution.
5. Orange G Solution staining is controlled under the microscope, subject to the removal of Acid Ponceau Solution. If used for samples containing active metal implants such as zinc magnesium, this step can be omitted, and then the bright green staining time can be extended to 10 min.