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Cat: G3661

Size: 100mL

Storage: RT, avoid light, valid for at least 1 year.

Introduction

Toluidine blue is a common synthetic dye, which belongs to quinone imine dyes. It has metachromatic property, so there are three kinds of dyeing solution configuration: acid, weak acid and alkaline. Weak acidity this method has a more balanced expression of heterochromatic coloring and natural coloring. It can be used to stain Nissl bodies and neurons at the same time. It can be used in paraffin section for general morphological analysis. It can also be used to differentiate islet cells.

Toluidine Blue O Solution, 1% in PBS can ensure the contrast coloring of background and heterochromatin at the same time due to the weak acidity of phosphate buffer.

Protocol (for reference only)

1. Fresh tissues are fixed with neutral formalin or paraformaldehyde, and routine dehydration embedding is performed.
2. Cut paraffin sections in 5 μ m thickness, then dewax paraffin sections and rehydrate in graded alcohol.
3. Preheat Toluidine Blue O Solution, 1% in PBS in a 50-60 °C incubator.
4. Stain with Toluidine Blue O Solution, 1% in PBS of 50-60°C for 20-40 min. Slightly wash with distilled water.
5. 95% ethanol differentiated rapidly for 3-5s. Absolute ethanol dehydration, transparent by xylene, seal with resinene.

Note

1. When using this reagent for the first time, it is recommended to take 1-2 samples for pre experiment.
2. When it was used to stain islet cells, the staining results were affected by the fixative. If Bouin fixative was used, only islet A1 cells were stained purple red. If conventional formalin fixative was used, both islet A1 and A2 cells were stained purple red.
3. For your safety and health, please wear experimental clothes and disposable gloves.

Result

Nissl Body	Purplish Blue
Nucleus	Light Blue
Background	Colorless or Light Blue

Related Products

- G3660 Toluidine blue staining Solution, For Cell
- G3662 Toluidine Blue O solution, 0.5% in PBS
- G3663 Toluidine Blue O Solution, 1% in Sodium Borate
- G3665 Toluidine Blue O Solution, 0.5% in Sodium Borate
- G3668 Toluidine Blue O Solution
- G3670 Mast Cells Stain Solution (Toluidine blue Method)