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**Cat:** G3710

**Size:** 2×50mL

**Storage:** 2-8°C, avoid light, valid for 6 months.

## Kit Components

Reagent	2×50mL	Storage
Reagent(A): Digestion Solution	50mL	2-8°C, avoid light
Reagent(B): Alcian Stain Solution	50mL	2-8°C, avoid light
Reagent(C): Negative Solution	10mL	RT

## Introduction

Alcian is a kind of copper titanium cyan conjugated dye, which was originally used for textile fiber dyeing. This kind of cationic dye can combine with acidic groups, that is, alcian blue forms an insoluble blue complex with anionic groups such as carboxyl and sulfate.

Hyaluronidase combined with alcian blue staining is used in the Hyaluronic Acid Stain Kit. The principle is that hyaluronidase can cleave the glycosidic bond of hyaluronic acid and chondroitin sulfate. Before alcian blue staining, the sample is pretreated with hyaluronidase. Compared with the untreated section, if the treated section is not stained, it indicates that there are hyaluronic acid and chondroitin sulfate in the tissue. If there is no effect, the tissue does not contain hyaluronic acid and chondroitin sulfate.

## Self Provided Materials

Series ethanol, Distilled Water, Warm Box

## Protocol(for reference only)

1. Divide two positive sections and two experimental sections two groups(Digestion Group and Control Group), dewax by xylene and series ethanol to water.
2. Treat Digestion Group with Digestion Solution and incubate at 37 °C for 3 h. Treat Control Group with Negative Solution and incubated at 37 °C for 3 h. Rinse with running water for 5 min.
3. Add alcian blue staining solution for 30 min, rinse with running water for 5 min.
4. Gradient ethanol dehydration, xylene transparent, mixed sealing agent sealing. Conventional fixation, conventional paraffin embedding.

## Result

Hyaluronic acid(Digestion Group)	Colorless
Hyaluronic acid(Control Group)	Blue

## Note

1. Positive sections should be used in each batch to detect enzyme activity.
2. There was a certain correlation between enzyme treatment efficiency and temperature, and 60 °C for 1 h could be used instead of 37 °C for 3 h.
3. For your safety and health, please wear lab clothes and disposable gloves.