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## Heidenhain's Iron-Hematoxylin Stain Kit

**Cat:** G4480

**Size:** 4×100mL

**Storage:** RT, avoid light, valid for 2 years.

### Kit Components

Reagent	4×100mL	Storage
Reagent(A): Heidenhain Differentiation	2×100mL	RT, avoid light
Reagent(B): Iron-Hematoxylin Solution	100mL	RT, avoid light
Reagent(C): Eosin Solution (Optional)	100mL	RT, avoid light

### Introduction

Hematoxylin-Eosin staining, which also named HE stain, is the most basic method of routine section staining in pathology. Hematoxylin is an alkaline natural dye, which can stain the nucleus. The main component of chromatin in the nucleus is DNA. In the double helix structure of DNA, the phosphate groups on the two nucleotide chains are outward, making the outer side of the double helix of DNA negatively charged and acidic. It is easy to dye with positively charged hematoxylin basic dye by ion bond or hydrogen bond.

Heidenhain's Iron-Hematoxylin Stain Kit use ammonium ferric sulfate as oxidant and differentiator, and can display different structures according to different differentiation degrees. After dyeing, all components are black or dark gray black. The hematoxylin staining of different tissue structure can be removed by Heidenhain differentiation solution at different speed, and the order of black color removal is: mitochondria, striated muscle, nuclear chromatin.

### Protocol (for reference only)

- 1) Fix tissue blocks, embed in paraffin and section.
- 2) Dewaxing and hydration of paraffin slices:
  - ① Dewax in xylene twice for 10mins each, Absolute ethanol twice for 2mins each.
  - ② 95% ethanol for 2mins, 85% ethanol for 2mins, 75% ethanol for 2mins, Distilled water for 2mins.
- 3) Mordantly stain with Heidenhain Differentiation for 1h (Cover dip dyeing), rinse with distilled water for 5-10s. (Refer Note 1)
- 4) Stain with Iron-Hematoxylin Solution for 1h (Suggested impregnation), Wash with tap water for 20-30s.
- 5) (Optional) Re-dyeing with Eosin Solution for 2-4min.
- 6) Differentiate with Heidenhain Differentiation or 1:1 water diluent of Heidenhain Differentiation. Wash alternately with tap water, the degree of differentiation was observed under microscope. (Refer Note 2)
- 7) Dehydration, transparency and sealing.
  - ① 75% ethanol for 10-20s, 85% ethanol for 10-20s
  - ② 95% ethanol twice for 1min, Absolute ethanol twice for 1min.
  - ③ Xylene twice for 1min. Seal with resin and view under the microscope.

### Result

Mitochondria, Striated Muscle, Myelin, Chromatin, etc	Dark Gray
---	-----------

### Note

1. The time of mordant dyeing by Heidenhain Differentiation and the time of staining by Iron-Hematoxylin Solution according to different fixative. Generally, the mordant dyeing and dyeing time can be controlled within 1h. For different fixatives, the reference time is: formalin, Bouin fixative, Carnoy fixative for 1h. Helly, Zenker and other dichromate fixative for 3h. osmium tetroxide, Flemming fixative for 24h.
2. Microscopically, the degree of differentiation is controlled until the desired structure appears. If the differentiation is excessive, hematoxylin can be used to re-dye and re-differentiate at the same time. In order to better control the degree of differentiation, we can also dilute Heidenhain differentiation with distilled water 2:1 before differentiation.
3. All the differentiation fluid should be washed thoroughly after the slice differentiation, otherwise the tissue is easy to fade. Slice dewaxing should be as clean as possible. Series ethanol should be replaced frequently, To prevent over staining, the dyeing time of frozen section must be short.
4. Re-dyeing the cytoplasmic (by Eosin or Orange G) can highlight nuclear chromatin, especially in the display of chromosomes or mitosis.
5. For your safety and health, please wear experimental clothes and disposable gloves.