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Pollen Viability Stain Solution (TTC Method)

Cat: G4800

Size: 10mL/50mL

Storage: 2-8°C, avoid light, valid for 1 year.

Introduction

The activity of pollen directly affects the process of pollination and fertilization, which is closely related to the yield of plants. Through the determination of pollen vigor, we can understand the development of pollen and master the morphological and physiological characteristics of sterile pollen. TTC is a redox pigment with a standard oxidation potential of 80 mV. A colorless solution is formed by dissolving in water, and a red insoluble triphenylmethane (TTF) is formed after reduction. TTC is widely used as hydrogen acceptor in enzyme experiments because it is stable and difficult to be oxidized. TTC reduction can express the activity of dehydrogenase, and then judge the activity of plant root or pollen.

Self Provided Materials

Slide, Incubator or Water Bath, Microscope

Protocol (for reference only)

1. Take fresh flowers that are about to open when they are ripe. Carefully remove petals and pistils to release pollen substances.
2. Put the pollen substances on the center of the slide, add 1-2 drops of Pollen Viability Stain Solution (TTC Method), cover the cover glass or take 0.5-1mL of Pollen Viability Stain Solution (TTC Method) and add it into a 1.5ml centrifuge tube. Add pollen material to the centrifuge tube until it is completely immersed in the staining solution and thoroughly mixed. (See Note 2)
3. Incubate at 37°C incubator for 30-45min. Observe under low power microscope, taking 5 visual fields from each slice. (See Note 3)

Result

Strong Vitality	Red
Low Vitality	Light Red
No Vitality or Sterility	Colourless

Calculation

Observe and count 100 pollens and calculate the percentage of active pollens.

The formula is: Percentage of pollen activity (%) = Number of active pollen / 100 × 100%

Note

1. After dyeing, observe under microscope immediately.
2. It is necessary to immerse the pollen completely in the staining solution.
3. The effect of 37°C constant temperature staining is significantly better than that of room temperature staining. If the staining is found to be light, the staining time can be extended to 1 hour.
4. For your safety and health, please wear experimental clothes and disposable gloves.