

## 50bp DNA Ladder instruction manual

**Item number:** M1800

**Specification:** 50T (250 $\mu$ L) /100T (500 $\mu$ L)

**Storage:** 2-8 $^{\circ}$ C valid for 6 months, -20 $^{\circ}$ C valid for 1 year.

### Product Introduction:

This product is composed of 8 strips of double-stranded DNA bands, suitable for the analysis of DNA bands in agarose gel electrophoresis.

This ready-to-use product contains 1 $\times$ loading buffer and takes 5 $\mu$ L directly for electrophoresis. It is easy to use and has clear electrophoretic image.

8 strips are divided into 50,100,150,200,250,300,400,500 bp, of which 250bp is 20ng/ $\mu$ L, and the rest is 10ng/ $\mu$ L.

### Composition of storage solution:

10mM Tris-HCl (pH8.4)


10mM EDTA

0.02% Bromophenol blue

5% glycerol

### Usage (for reference only) :

1. Take 5 $\mu$ L of this product and add it into the sample hole of agar-agar gel (add 1 $\mu$ L for each 1mm sample hole width, if the sample hole is wider, the sample amount can be appropriately increased) for electrophoresis.
2. It is recommended that the gel concentration is 2-3% agarose gel, the electrophoretic voltage is 4-10v/cm, and the electrophoretic time is 30-40min.
3. Nucleic acid stain and observe electrophoretic bands under UV lamp.



5 $\mu$ l loading, 2.5% agarose gel electrophoresis diagram

### Precautions:

1. Change the electrophoresis buffer and use the newly formulated agarose gel in time to avoid affecting the electrophoresis result.
2. Since nucleic acid binding dyes can affect the migration of DNA during electrophoresis, it is recommended to perform gel bliscer after electrophoresis.

**Related products:**

<i>A8201</i>	<i>agarose</i>
<i>D1020</i>	<i>10 x DNA loading buffer</i>
<i>G8142</i>	<i>GoldView Type II nucleic acid stain (5000×)</i>
<i>T1060</i>	<i>50 x TAE buffer</i>
<i>T1050</i>	<i>5 x TBE buffer</i>

**Note:** For more information about this product, please refer to the Solarbio website.