

PCR MasterMix instruction manual

Item	No.Name	Specification
PC1120	2 x Taq PCR MasterMix (with Dye)	1mL /1mL×5
PC1150	2 x Taq PCR MasterMix (without dye)	1mL /1mL×5
PC1320	2×Pfu PCR MasterMix (including dye)	1mL /1mL×5
PC1350	2×Pfu PCR MasterMix (without dye)	1mL /1mL×5
PC1220	2×Taq Plus MasterMix (with dye)	1mL /1mL×5
PC1250	2×Taq Plus MasterMix (without dye)	1mL /1mL×5

Storage: Store at -20°C for 12 months, store at 2-8°C for a short period of one week to avoid repeated freezing and thawing.

Product introduction:

The PCR MasterMix produced by Solarbio company contains three enzymes (Taq, Pfu, Taq plus), each enzyme corresponds to two kinds of PCR MasterMix (dye and dye free), a total of six kinds of PCR MasterMix.

Taq: The highest amplification efficiency of heat-resistant DNA polymerase, its amplification rate of about 1-2kb/min. The amplification base error rate is about 10^{-5} . The product has A on the end and can be directly used for TA cloning.

Pfu: Currently the highest fidelity of heat-resistant DNA polymerase, base error rate is 10^{-6} , but the amplification efficiency is lower than Taq enzyme, generally can very well amplify the fragment below 2kb. Its amplification speed is about 500bp/min. The product is A flat end and can be used for TA cloning only after adding A.

Taq Plus: An equal mixture of Taq and Pfu. A combination of high amplification efficiency and high fidelity. Amplification efficiency is higher than Pfu, fidelity is better than Taq. The amplification rate was about 1000bp/min. The product with A on the end can be directly used for TA cloning.

After the reaction of dye-containing PCR MasterMix, it can be directly detected by electrophoresis. After the reaction of dye-free PCR MasterMix, the sample buffer is added for electrophoretic detection.

Product composition (2×) :

20mM	Tris-HCl (pH 8.3)
100mM	KCl
3mM	MgCl ₂
500μM	dNTP each
0.1U	Polymerase/μL
ddH ₂ O, other stabilizers and reinforcers	

Application examples (for reference only) :

1. Prepare 50μL PCR reaction system

Template <	1μg
primer1(10μM)	2μL
Primer2(10μM)	2μL
2×MasterMix	25μL
ddH ₂ O	up to 50μL

2. PCR reaction cycle setting

94°C	3min	} 30 cycles
94°C	30s	
55°C	30s	
72°C	500-2000bp/1min	
72°C	5min	

3 Result detection: After the reaction, take 5µL reaction product, agarose gel electrophoresis detection.

Related products:

<i>A8201</i>	<i>Agarose</i>
<i>T1060</i>	<i>50 x TAE buffer</i>
<i>D1010</i>	<i>6×DNA Loading Buffer</i>
<i>M1070</i>	<i>D2000 plus DNA Ladder</i>
<i>PC1100</i>	<i>Taq DNA Polymerase</i>
<i>PC2100</i>	<i>dNTPs Mix(2.5mM each)</i>
<i>PC2200</i>	<i>dNTPs Mix(10mM each)</i>
<i>RP1100</i>	<i>General Purpose RT-PCR Kit (M-MLV)</i>
<i>SR1110</i>	<i>2×SYBR Green PCR Mastermix</i>
<i>SR2110</i>	<i>2×Taqman PCR MasterMix</i>