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Cat: G1136

Size: 2×100mL

Storage: 2-8°C, avoid light, valid for 6 months.

Kit Components

Reagent	2×100mL	Storage	
Reagent(A): Tannic Acid Working Solution	A1: Tannic Acid Diluent	50mL	RT
	A2: Tannic Acid Solution	50mL	RT, avoid light
Mix A1 and A2 with the ratio of 1:1 to form Tannic Acid Working Solution before use.			
Reagent(B): Ammonia Silver Staining Solution	100mL	2-8°C, avoid light	

Introduction

Bacterial flagellum is the motor organ of bacteria. *Helicobacter pylori* can pass through the mucus layer of gastric epithelial cells from the strongly acidic gastric cavity to the neutral environment of gastric epithelial cells, which is a good example of flagellum movement. Flagella morphology, quantity and location of flagella in bacteria can be observed by flagella staining. Flagellum quantity and its distribution on bacteria are one of the important bases for identifying bacteria.

This kit uses silver staining method. The advantage of the method is that Ammonia Silver Staining Solution is used as the core dye. The reagent is more sensitive, the operation is simple and the result is more reliable.

Self Provided Materials

Alcohol lamp, slide, distilled water, inoculation ring, microscope.

Protocol (for reference only)

- Drop 2 drops of distilled water on a clean, fat-free slide.
- Select the sterile distilled water with the inoculation ring, then contact the colony on the blood plate, allow the bacteria to swim into the inoculation ring distilled water, and then move the inoculation ring to the top of the distilled water on the slide twice.
- Shake the slide to distribute the bacteria. Do not grind and stir with the slide to avoid flagella falling off.
- Dry in a room temperature or 37°C incubator (not fixed by flame).
- Add the Tannic Acid Working Solution to a clean container or test tube, mix it well, heat it slowly with an alcohol lamp, cool it slightly, and dye it immediately.
- Drop the heated Tannic Acid Working Solution onto the slide and dye it for 30-40 s. Rinse it slowly with distilled water and shake it dry.
- Drop Ammonia Silver Staining Solution, heat in small flame until bubbling, rinse with distilled water, and dry naturally.
- Microscopic examination should begin at the edge of the smear and gradually move to the center from the outside to the inside. Where bacteria are less distributed, flagellum is easy to observe. Where bacteria are concentrated, flagellum is blocked by fungi and is not easy to observe.

Result

Thalli	Dark Brown
Flagellum	Brown

Note

- Slides should be clean and free from oil pollution.
- The bacteria sample for staining must be in the active stage of growth, which should be successive subculture for many times.
- Operate carefully in the dyeing process to prevent flagella falling off. It is not appropriate to fix with high-temperature flame.
- The prepared Tannic Acid Working Solution should be used within 10 min as far as possible.
- Ammonia Silver Staining Solution is unstable. It should be stored avoiding light at 4°C strictly and should be discarded when serious turbidity occurs.
- For your safety and health, please wear experimental clothes and disposable gloves.