

# Plant Soluble Sugar Content Assay Kit

**Note:** Take two or three different samples for prediction before test.

**Detection instrument:** Spectrophotometer/Microplate Reader

**Cat No:** BC0035

**Size:** 100T/96S

## Components:

Reagent I: Powder×2. Storage at 4°C and protect from light.

Reagent II: Liquid 5 mL×1. Storage at 4°C.

Standard: Powder×1, containing 10 mg of anhydrous glucose(loss on drying<0.2%). Storage at 4°C. Dissolve the standard with 1 mL of distilled water to generate a 10 mg/mL glucose standard solution, store at 4°C and use within one week. Or dissolve the standard with saturated benzoic acid solution for a longer time.

Preparation of working solution: Add 2.5 mL of Reagent II to Reagent I, fully dissolved before use. If it's hard to dissolve, it can be heated and stirred to auxiliary dissolved. The reagent that can not be used up can be kept at 4 °C for one week.

## Product Description:

Sugar is one of the important components of plant, also is the main material of metabolism and storage. Total sugar refers to the reductive monosaccharide in the sample and the sucrose, maltose and starch which can be hydrolyzed or partially hydrolyzed to glucose by this test method.

Anthrone colorimetric method can be used for the determination of soluble monosaccharides, oligosaccharides and polysaccharides. It has the advantages of high sensitivity, simple and quick, and suitable for the determination of trace samples.

## Reagents and Equipment Required but Not Provided:

Spectrophotometer/microplate reader, centrifuge, micro glass cuvette/96 well flat-bottom plate, transferpettor, concentrated sulfuric acid, mortar/homogenizer, distilled water.

## Procedure:

### I. Sample Extraction:

Suggested 0.1-0.2 g of sample with 1 mL of distilled water, fully grinding, boiling water bath for 10 minutes (cover tightly to prevent water loss). After cooling, centrifuge at 8000×g at room temperature for 10 minutes, then take the supernatant in 10 mL tube, add distilled water to 10 mL, mix thoroughly.

### II. Determination procedure:

1. Preheat the spectrophotometer/microplate reader 30 minutes, adjust wavelength to 620 nm, set zero with distilled water.
2. Set the temperature of water bath to 95°C.
3. Standard working solution: Dilute the 10 mg/mL glucose standard solution to 0.3, 0.2, 0.1, 0.05, 0.025,

0.0125 mg/mL with distilled water.

4. Add reagents with the following list(reaction in EP tube):

Reagent (μL)	Blank tube (B)	Test tube (T)	Standard tube (S)
Sample	-	40	-
Standard	-	-	40
Distilled water	80	40	40
Working solution	20	20	20
Concentrated sulfuric acid	200	200	200

Mix thoroughly, put the reaction solution in 95°C water bath for 10 minutes (cover tightly to prevent water loss), detect the absorbance after cooling to room temperature. Transfer 200 μL to a micro cuvette or 96 well flat-bottom plate. Measure the absorption value at 620 nm.  $A_{\text{blank tube}}=A(B)$ ,  $A_{\text{test tube}}=A(T)$ ,  $A_{\text{Standard tube}}=A(S)$ .  $\Delta A = A(T)-A(B)$ ,  $A = A(S)-A(B)$ . The blank tube is needed test once or twice.

### III. Calculation:

1. According to the standard curve, add  $\Delta A$  into the formula (x) to calculate the sample concentration y (mg/mL).

2. Sample weight

$$\text{Soluble sugar(mg/g fresh weight)}=y \times V1 \div (W \times V1 \div V2)=10 \times y \div W$$

3. Protein concentration

$$\text{Soluble sugar(mg/mg prot)}=y \times V1 \div (V1 \times Cpr)=y \div Cpr$$

V1: Sample volume, 0.04 mL;

V2: Extraction volume, 10 mL;

Cpr: Sample protein concentration, mg/mL;

W: Sample weight, g.

### Note:

1. Dilute the sample and measure again if  $\Delta A$  is greater than light absorption value in linear range. When calculation, multiply the calculation formula by the corresponding dilution factor.

2. Due to the strong corrosiveness of concentrated sulfuric acid, please opera with caution.

### Recent Product citations:

[1] Li B, Ding Y, Tang X, et al. Effect of L-Arginine on Maintaining Storage Quality of the White Button Mushroom (*Agaricus bisporus*)[J]. Food and Bioprocess Technology, 2019, 12(4): 563-574.

[2] Zeyong Zhang, Huanhuan Liu, Ce Sun, et al. A C<sub>2</sub>H<sub>2</sub> zinc-finger protein OsZFP213 interacts with OsMAPK3 to enhance salt tolerance in rice. Journal of Plant Physiology. October 2018;(IF2.825)

### References:

[1] Buysse J A N, Merckx R. An improved colorimetric method to quantify sugar content of plant tissue[J]. Journal of Experimental Botany, 1993, 44(10): 1627-1629.

[2] Bodelón O G, Blanch M, Sanchez-Ballesta M T, et al. The effects of high CO<sub>2</sub> levels on anthocyanin composition, antioxidant

**Related products:**

BC0330/BC0335 Trehalose Assay Kit

BC0340/BC0345 Glycogen Assay Kit

BC2530/BC2535 Sorbitol Dehydrogenase Activity Assay Kit

BC0230/BC0235 Reducing Sugar Content Assay Kit

**Technical Specifications:**

Minimum Detection Limit: 0.0025 mg/mL

Linear Range: 0.003125-0.4 mg/mL