

Ferrous Ion Content Assay Kit

Note: It is necessary to predict 2-3 large difference samples before the formal determination.

Operation Equipment: Spectrophotometer/microplate reader

Cat No: BC5415

Size:100T/96S

Components:

Reagent I: Liquid 125 mL×1. Store at 2-8°C.

Reagent II: Liquid 13 mL×1. Store at 2-8°C.

Standard: Powder×1, 10 mg $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$. Store at 2-8°C. Add 900 μL of distilled water and 20 μL of concentrated sulfuric acid before use and shake to dissolve. The Fe^{2+} standard solution of 40 mmol/L could be stored at 2-8°C for two weeks.

Product Description:

Iron is one of the essential trace elements in human body. Fe^{2+} is the main component of hemoglobin, myoglobin, cytochrome and other enzyme systems, which could assist in the transport of oxygen and promote fat oxidation. Iron deficiency can easily cause anemia, metabolic disorders, and affect the immune function of the body.

Fe^{2+} could react with Tripyridyltriazine to form a kind of blue compound under acid condition, which has an absorption peak at 593 nm. Changes of the absorbance at 593nm could be measured to reflect ferrous ion content.

Reagents and Equipment Required but Not Provided:

Spectrophotometer/microplate reader, balance, low temperature table centrifuge, constant temperature incubator/water bath, micro glass cuvette/96 well flat-bottom plate, pipette, mortar/homogenizer/cell ultrasonic crusher, ice, distilled water and **concentrated sulfuric acid**.

Procedure:

I. Sample preparation:

- Tissue:** according to tissue weight (g): Reagent I volume (mL) is 1:5-10. (It is recommended that add 1 mL of Reagent I to 0.1 g tissue). Homogenate in ice bath, then centrifuge at 10000 g for 10 minutes at 4°C. Take the supernatant for test.
- Bacteria/cells:** according to the number of bacteria/cells (10^4): the volume of Reagent I (mL) is 500~1000:1. It is suggested that add 1 mL of Reagent I to 500 million of cells. Breaking bacteria/cells by ultrasonic wave in ice bath (power 200W, ultrasonic 3s, interval 7s, total time 5 min). Centrifuge at 10000 g 4°C for 10 minutes. Take the supernatant on ice for test.
- Serum (plasma) or other liquid samples:** detect directly. Centrifuge before detecting if there are precipitation in the liquid.

II. Determination procedure:

- Preheat spectrophotometer/microplate reader for 30 minutes, adjust wavelength to 593 nm, set spectrophotometer counter to zero with distilled water.

- Standard working solution: Prepare 400 $\mu\text{mol/L}$ standard solution with 10 μL of 40 mmol/L standard solution and 990 μL of distilled water. Dilute 400 $\mu\text{mol/L}$ standard solution with Reagent I to 50, 25, 12.5, 6.25, 3.125, 1.5625, 0.78125 $\mu\text{mol/L}$ for standby.
- Add reagents with the following list:

Reagent (μL)	Test tube (A_T)	Standard tube (A_S)	Blank tube (A_B)
Sample	200	-	-
Standard solution	-	200	-
Reagent I	-	-	200
Reagent II	100	100	100
Mix thoroughly. React at 37°C for 10 minutes.			
Trichloromethane	100	-	-
Mix thoroughly for 5 minutes and centrifuge at 12000 g for 10 minutes at room temperature. Take 200 μL upper inorganic phase solution in micro glass cuvette/96 well flat-bottom plate. Measure absorbance at 593 nm, recorded as A_T . $\Delta A_T = A_T - A_B$.		Measure absorbance at 593 nm, recorded as A_B , and A_S . $\Delta A_S = A_S - A_B$. Blank tube and standard curve only need to test once or twice.	

III. Ferrous Ion Content Calculations

1. Standard curve

Taking the concentration of each standard solution as the x-axis and its corresponding ΔA_S as the y-axis, draw a standard curve to get the standard equation $y = kx + b$, and bring ΔA_T into the equation to get x ($\mu\text{mol/L}$).

2. Calculation

- Serum (plasma) or other liquid samples: Ferrous ion content ($\mu\text{mol/L}$) = x
- Protein concentration: Ferrous ion content ($\mu\text{mol/mg prot}$) = $x \times 10^{-3} \times V_E \div (C_{pr} \times V_E) = 0.001x \div C_{pr}$
- Sample weight: Ferrous ion content ($\mu\text{mol/g weight}$) = $x \times 10^{-3} \times V_E \div W = 0.001x \div W$
- Bacteria/cells number: Ferrous ion content ($\mu\text{mol}/10^6 \text{ cell}$) = $x \times 10^{-3} \times V_E \div N = 0.001x \div N$

C_{pr} : Supernatant sample protein concentration, mg/mL ;

V_E : Added Reagent I volume, 1 mL;

W: Sample weight, g;

N: Total number of bacteria or cells, per 10^6 ;

10^{-3} : Unit conversion factor, 1 $\mu\text{mol/L} = 10^{-3} \mu\text{mol/mL}$.

Note:

- It is better to prepare standard solution before using because standard solution diluted with Reagent I easily fail.
- If A_T is close to A_B or ΔA_T is too low, it is recommended to increase sample supernatant size before determination. If $A_T > 0.4$, it is recommended to dilute sample supernatant with Reagent I before determination. And modify the calculation formula.

3. It is suggested to take upper inorganic phase solution carefully to avoid taking lower trichloromethane solution.

Experimental example:

1. Take 0.1055g mouse liver, add 1 mL of Reagent I, grind the homogenate with ice bath. Then operate according to the determination steps, calculate $\Delta A_T = A_T - A_B = 0.272 - 0.050 = 0.222$. Bring the result into the standard curve $y = 0.0076x - 0.0013$ and calculate $x = 29.38$. The result is calculated according to sample weight:

Ferrous ion content ($\mu\text{mol/g weight}$) $= 0.001x \div W = 0.029 \mu\text{mol/g weight}$.

2. Take 200 μL of calf serum, operate according to the determination steps, calculate $\Delta A_T = A_T - A_B = 0.301 - 0.050 = 0.251$. Bring the result into the standard curve $y = 0.0076x - 0.0013$ and calculate $x = 33.20$. The result is calculated according to liquid volume:

Ferrous ion content ($\mu\text{mol/L}$) $= x = 33.20 \mu\text{mol/L}$.

References:

[1] Dong C, Yang M, Wang W. Study on spectrophotometric determination of Fe(II) and Fe(III) with 2, 4, 6- tri(2'-pyridyl)-1, 3, 5-triazine. [J]. Chinese Journal of Analysis Laboratory, 2004, (01):76-78.

[2] Huang Y, Ma H, Xu J, et al. Development and Validation of Reference Methods for Determination of Serum iron. [J]. Chinese Journal of Laboratory Diagnosis, 2011, 15(03):453-457.

[3] Wang H, Liu B, Ding Z, et al. Ferene method flow injection analysis as optimized in situ analysis of dissolved iron in marine waters. [J]. Marine Sciences, 2016, 40(05):82-87.

Related Products :

BC1730/BC1735 Serum Ferri Ion Content Assay Kit

BC2810/BC2815 Blood Zinc Content Assay Kit

BC2860/BC2865 Serum Total Iron Binding Capacity(TIBC) Assay Kit

BC4350/BC4355 Tissue Iron Content Assay Kit

