

50×Denhardt's Solution

Cat: D1080

Size: 50mL

Storage: Store at -20°C, valid for 1 year.

Introduction:

Solarbio produced 50×Denhardt's solution is mainly composed of Ficoll 400, PVP, BSA, mainly used for DNA hybridization of single copy sequence, DNA fixed on nylon membrane hybridization. Solarbio's Denhardt's solution has been treated aseptically.

Instructions for use:

Follow the specific requirements of the experiment, or dilute in hybridization buffer(usually SSC containing SDS, denatured salmon extract DNA)

Notes:

1. If the amount used each time is small, it can be used after appropriate packaging to avoid repeated freezing and thawing.
2. Pay attention to aseptic operation to prevent bacterial contamination.
3. For your safety and health, please wear a lab coat and gloves.

Related Products:

H1060 Salmon Essence DNA, 10mg/mL

F8120 Formamide

S1035 20×SSC, PH5.3

S1090 20×SSC, PH7.4

S1160 20×SSPE, PH7.4

D2810 DNA Annealing Buffer(5×)