


2019-nCoV S2 Protein

Catalog#:P00308 Derived from E.coli

DESCRIPTION	2019-nCoV S2 Protein is produced by our E.coli expression system with His tag. Uniprot: P0DTC2 Gene ID: 43740568
Size	32 kDa
FORMULATION	Lyophilized from a 0.2 um filtered solution of PBS, PH 7.4.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 95% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding Lys301-Thr588 is expressed with a 6*His tag at the C-terminus.
BACKGROUND	/



SARS-CoV-2(2019-nCoV))Spike S2-His Recombinant Protein(2)