


Recombin human CA2 protein

Catalog#:P02392 Derived from E.coli

DESCRIPTION	Recombinant Human CA2 protein is produced by our E.coli expression system with His tag. Uniprot: P00918 Gene ID: 760
Size	29 KDa
FORMULATION	Lyophilized from a 0.2 µm filtered solution of Tris-HCl, PH 8.0.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 85% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding 1-138aa is expressed with a 6*His tag at the N-terminus.
BACKGROUND	The protein encoded by this gene is one of several isozymes of carbonic anhydrase, which catalyzes reversible hydration of carbon dioxide. Defects in this enzyme are associated with osteopetrosis and renal tubular acidosis. Two transcript variants encoding different isoforms have been found for this gene.



Recombin human CA2 protein