

Recombin human CAT protein

Catalog#:P02463 Derived from E.coli

DESCRIPTION	Recombinant Human CAT protein is produced by our E.coli expression system with His tag. Uniprot: P04040 Gene ID: 847
Size	37 KDa
FORMULATION	Lyophilized from a 0.2 µm filtered solution of Tris-HCl, PH 8.0.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 85% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding 200-527aa is expressed with a 6*His tag at the N-terminus.
BACKGROUND	This gene encodes a nuclear pre-mRNA splicing factor. The encoded protein specifically recognizes the intron branch point sequence at the 3' splice site, together with the large subunit of U2 auxiliary factor (U2AF), and is required for the early stages of spliceosome assembly. It also plays a role in nuclear pre-mRNA retention and transcriptional repression. The encoded protein contains an N-terminal U2AF ligand motif, a central hnRNP K homology motif and quaking 2 region which bind a key branch-site adenosine within the branch point sequence, a zinc knuckles domain, and a C-terminal proline-rich domain. Alternative splicing results in multiple transcript variants.

