

Recombin human LPXN protein

Catalog#:P02897 Derived from E.coli

DESCRIPTION	Recombinant Human LPXN protein is produced by our E.coli expression system with His tag. Uniprot: O60711 Gene ID: 9404
Size	43 KDa
FORMULATION	Lyophilized from a 0.2 µm filtered solution of Tris-HCl, PH 8.0.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 85% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding 1-386aa is expressed with a 6*His tag at the N-terminus.
BACKGROUND	The product encoded by this gene is preferentially expressed in hematopoietic cells and belongs to the paxillin protein family. Similar to other members of this focal-adhesion-associated adaptor-protein family, it has four leucine-rich LD-motifs in the N-terminus and four LIM domains in the C-terminus. It may function in cell type-specific signaling by associating with PYK2, a member of focal adhesion kinase family. As a substrate for a tyrosine kinase in lymphoid cells, this protein may also function in, and be regulated by, tyrosine kinase activity. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

