


Recombinant human CAV1 protein

Catalog#:P04501 Derived from E.coli

DESCRIPTION	Recombinant Human CAV1 protein is produced by our E.coli expression system with His tag. Uniprot: Q03135 Gene ID: 857
Size	20 kDa
FORMULATION	Lyophilized from a 0.2 µm filtered solution of Tris-HCl, PH 8.0.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 85% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding 1-178aa is expressed with a 6*His tag at the N-terminus.
BACKGROUND	The scaffolding protein encoded by this gene is the main component of the caveolae plasma membranes found in most cell types. The protein links integrin subunits to the tyrosine kinase FYN, an initiating step in coupling integrins to the Ras-ERK pathway and promoting cell cycle progression. The gene is a tumor suppressor gene candidate and a negative regulator of the Ras-p42/44 mitogen-activated kinase cascade. Caveolin 1 and caveolin 2 are located next to each other on chromosome 7 and express colocalizing proteins that form a stable hetero-oligomeric complex. Mutations in this gene have been associated with Berardinelli-Seip congenital lipodystrophy. Alternatively spliced transcripts encode alpha and beta isoforms of caveolin 1.



Recombinant human CAV1 protein