

Recombinant human BMP7 protein

Catalog#:P05985 Derived from E.coli

DESCRIPTION	Recombinant Human BMP7 protein is produced by our E.coli expression system with His tag. Uniprot: P18075 Gene ID: 655
Size	40 kDa
FORMULATION	Lyophilized from a 0.2 µm filtered solution of Tris-HCl, PH 8.0.
SHIPPING	The product is shipped at -20°C temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Reconstituted protein solution can be stored at 4-7°C for 2-3 months, stable at < -20°C for 1-2 years.
RECONSTITUTION	Reconstituted protein solution can be diluted with distilled PBS. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	(Liquid)Concentration: 1 mg/mL as determined by BCA. Purity: 85% as determined by reducing SDS-PAGE.
AMINOACID	Recombinant protein is produced by our E.coli expression system and the target gene encoding 51-400aa is expressed with a 6*His tag at the N-terminus.
BACKGROUND	This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which plays a role in bone, kidney and brown adipose tissue development. Additionally, this protein induces ectopic bone formation and may promote fracture healing in human patients.

