

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

## **Recombinant Mouse RANK L**

Catalog#:P02067 Derived from Human Cells

	Recombinant Mouse TNF-related Activation-induced Cytokine is produced by our Mammalian expression system and the target gene encoding Arg72-Asp316 is expressed with a 6His tag at the N-terminus.
DESCRIPTION	Accession#: BAA97257.1
	<b>Known as:</b> Tumor necrosis factor ligand superfamily member 11; Tnfsf11;
	Osteoclast differentiation factor; ODF; Osteoprotegerin ligand; OPGL; Receptor
	activator of nuclear factor kappa-B ligand; RANKL; TNF-related
	activation-induced cytokine; TRANCE; CD254
FORMULATION	Lyophilized from a 0.2µm filtered solution of PBS, pH 7.4.
SHIPPING	The product is shipped at ambient temperature.
	Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at <-20°C, though stable at room
	temperature for 3 weeks.
	Reconstituted protein solution can be stored at 4-7°C for 2-7 days.
	Aliquots of reconstituted samples are stable at < -20°C for 3 months.
	Always centrifuge tubes before opening. Do not mix by vortex or pipetting.
RECONSTITUTION	It is not recommended to reconstitute to a concentration less than 100µg/ml.
RECONSTITUTION	Dissolve the lyophilized protein in distilled water.
	Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY	Mol Mass:28.3kDa AP Mol Mass:30-35kDa, reducing conditions.
	<b>Purity</b> : Greater than 95% as determined by reducing SDS-PAGE.
CONTROL	Endotoxin: Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
	Mouse tumor necrosis factor ligand superfamily member 11(Tnfsf11) is a member of the tumor necrosis factor (TNF) cytokine family. Tnfsf11 is widely expressed in cells including T cells and T cell rich organs, such as the truly conditions of the tumor necrosis factor (TNF) cytokine family. Tnfsf11 is widely expressed in cells including T cells and T cell rich organs, such as the truly conditions of the tumor necrosis factor ligand superfamily member 11(Tnfsf11) is a member of the tumor necrosis factor (TNF) cytokine family. Tnfsf11 is widely expressed in cells including T cells and T cell rich organs, such as the tumor necrosis factor (TNF) cytokine family.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	and TNFRSF11A/RANK. This cytokine can bind to TNFRSF11B/OPG and TNFRSF11A/RANK. This fill is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.    No.   MK   R   120   60   60   60   60   60   60   60
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.   KDB MK R  120 90 60 40  SDS-PAGE
BACKGROUND	lymph nodes. This cytokine can bind to TNFRSF11B/OPG andTNFRSF11A/RANK. Tnfsf11 is involved in a number of fundamental biological processes such as acting as regulator of interactions between T-cells and dendritic cells, the regulation of the T-cell-dependent immune response and enhancing bone-resorption in humoral hypercalcemia of malignancy. It augments the ability of dendritic cells to stimulate naive T-cell proliferation.   KDB MK R  120 90 60 40  SDS-PAGE