

Immune Serum

Summary:

Quantity size: 0.5mL/1mL/5mL

Immunogen: IgM

Double agar diffusion assay: 1:4-32

Western bloting: 1:500-10000

ELISA: 1:500-50000

Storage: Shipped at 4°C, Store at -20°C at least one year (Avoid repeated freeze/thaw cycles).

Background:

Antisera against whole serums are obtained by immunizing host animals with whole serum. This product is reddish transparent liquid, containing 0.03% Proclin300.

Related products:

SA201	Rabbit anti-Human IgG(Immune Serum)	A8201	Agarose
SA231	Rabbit anti-Mouse IgG(Immune Serum)	SP034	Rabbit IgG
SA235	Rabbit anti-Chicken IgG(Immune Serum)	SP031	Mouse IgG
SPA101	Goat anti-Human IgG(Purification Antibody)	SP037	Pig IgG
SPA131	Goat anti-Mouse IgG(Purification Antibody)	SE237	Rabbit Anti-pig IgG-HRP
SF134	Goat anti-Rabbit IgG-FITC	SE101	Goat anti-Human IgG-HRP
SR134	Goat anti-Rabbit IgG-RBITC	SHB131	Goat anti-Mouse IgG-Biotin
SHB134	Goat anti-Rabbit IgG-Biotin		

Related Articles:

[1] Lili Tan,Limin Liu,Zhenyu Jiang,et al. Inhibition of microRNA-17-5p reduces the inflammation and lipid accumula -tion, and up-regulates ATP-binding cassette transporterA1 in atherosclerosis. Journal of Pharmacological Sciences. April 2019;139(4):280-288. (IF 2.439)

Important Note: Potency is subject to the product label. Optimal working dilutions must be determined by the end user.