

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

Recombinant Human IL-13

Catalog#:P01686 Derived from Human Cells

	Catalog#:P01686 Derived from Human Cells		
DESCRIPTION	Recombinant Human Interleukin-13 is produced by our Mammalian expression		
	system and the target gene encoding Gly21 -Asn132 is expressed with a Fc tag at		
	the C-terminus.		
	Accession#: AAK53823		
	Known as: Interleukin-13; IL-13		
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.		
SHIPPING	The product is shipped at ambient temperature.		
	Upon receipt, store it immediately at the temperature listed below.		
STORAGE	Lyophilized protein should be stored at <-20°C, though stable at room		
	temperature for 3 weeks.		
	Reconstituted protein solution can be stored at 4-7°C for 2-7 days.		
	Aliquots of reconstituted samples are stable at < -20°C for 3 months.		
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting.		
	It is not recommended to reconstitute to a concentration less than 100µg/ml.		
	Dissolve the lyophilized protein in distilled water.		
	Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
QUALITY	Mol Mass:39.2kDa AP Mol Mass:39-60kDa, reducing conditions.		
CONTROL	Purity: Greater than 95% as determined by reducing SDS-PAGE.		
CONTROL	Endotoxin: Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.		
BACKGROUND	Interleukin-13 is also known as IL-13. It is a protein that in humans is encoded by the IL13 gene. Interleukin-13 is an immunoregulatory cytokine produced primarily by activated Th2 cells. It is involved in several stages of B-cell maturation and differentiation. It up-regulates CD23 and MHC class II expression, and promotes IgE isotype switchin of B cells. This cytokine down-regulates macrophage activity, thereby inhibits the production of pro-inflammatory cytokines and chemokines. This cytokine is found to be critical to the pathogenesis of allergen-induced asthma but operates through mechanisms independent of IgE and eosinophils.		
	kDa MK R		
120 90 60 40 SDS-PAGE 30			
			20