

Recombinant Human IL-4

Catalog#:P01499 Derived from Human Cells	
Recombinant Human Interleukin-4 is produced by our Mammalian expression	
DESCRIPTION	system and the target gene encoding His25-Ser153 is expressed with a Fc tag at
	the C-terminus.
	Accession#: P05112
	Known as: Interleukin-4; IL-4; B-Cell Stimulatory Factor 1; BSF-1; Binetrakin;
	Lymphocyte Stimulatory Factor 1; Pitrakinra; IL4
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
SHIPPING	The product is shipped at ambient temperature.
	Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at <-20°C, though stable at room
	temperature for 3 weeks.
	Reconstituted protein solution can be stored at 4-7°C for 2-7 days.
	Aliquots of reconstituted samples are stable at $< -20^{\circ}$ C for 3 months.
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting.
	It is not recommended to reconstitute to a concentration less than 100µg/ml.
	Dissolve the lyophilized protein in distilled water.
	Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
	Mol Mass:41.9kDa AP Mol Mass:40-55kDa, reducing conditions.
QUALITY	Purity : Greater than 95% as determined by reducing SDS-PAGE.
CONTROL	Endotoxin : Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL test.
BACKGROUND	Interleukin-4 (IL-4) is a pleiotropic cytokine that regulates diverse T and B cell responses including cell proliferation, survival and gene expression. IL-4 is produced by mast cells, T cells, and bone marrow stromal cells. IL-4 regulates the differentiation of naive CD4+ T cells into helper Th2 cells, characterized by their cytokine-secretion profile that includes secretion of IL-4, IL-5, IL-6, IL-10, and IL-13, which favor a humoral immune response. Another dominant function of IL-4 is the regulation of immunoglobulin class switching to the IgG1 and IgE isotypes. Excessive IL-4 production by Th2 cells has been associated with elevated IgE production and allergic response.
	kDa MK R
120 90	
	60
	40
	SDS DACE 30
	SDS-PAGE
	20
	14