

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

## **Recombinant Mouse SOST**

Catalog#:P02126 Derived from Human Cells

Catalogn. 102120 Derived from Human Cens	
DESCRIPTION	Recombinant Mouse Sclerostin is produced by our Mammalian expression system and the target gene encoding Gln24-Tyr211 is expressed with a 6His tag at the C-terminus.  Accession#: AAK13455.1  Known as: Sclerostin; Sost
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
SHIPPING	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at<-20°C, though stable at room temperature for 3 weeks.  Reconstituted protein solution can be stored at 4-7°C for 2-7 days.  Aliquots of reconstituted samples are stable at < -20°C for 3 months.
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY	Mol Mass:21.9kDa AP Mol Mass:32kDa, reducing conditions.  Purity: Greater than 95% as determined by reducing SDS-PAGE.
CONTROL	Endotoxin: Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
BACKGROUND	Sclerostin (SOST) is a member of the cerberus/DAN family, a group of secreted glycoproteins characterized by a cysteine-knot motif. Cerberus/DAN family members are putative BMP antagonists, and include Dan, Cerberus, Gremlin, PRDC, and Caronte. While the overall sequence identity between members of the family is low, they have conserved spacing of six cysteine residues. Cerberus and Dan have an additional cysteine residue used for dimerization; However, SOST does not and is secreted as a monomer. SOST was originally identified as an important regulator of bone homeostasis. SOST is expressed by osteoclasts in developing bones of mouse embryos, including both intramembranously forming skull bones and endochondrally forming long bones. SOST plays a physiological role as a negative regulator of bone formation by repressing BMP-induced osteogenesis. SOST has been shown to have unique ligand specificity, binding
	BMP-5, -6, and -/ with high affinity and BMP-2 and -4 with low affinity.
LDa MAY P	

