

Recombinant Human PSMA

Catalog#:P01843 Derived from Human Cells	
DESCRIPTION	 Recombinant Human Glutamate Carboxypeptidase 2 is produced by our Mammalian expression system and the target gene encoding Lys44-Ala750 is expressed with a 6His tag at the N-terminus. Accession#: Q04609 Known as: Glutamate carboxypeptidase 2; FGCP; GCPII; mGCP; NAALADase I; PSMA; Cell growth- inhibiting gene 27 protein; Folate hydrolase 1
FORMULATION	Supplied as a 0.2µm filtered solution of 20mM MES, 150mM NaCl, 5% Trehalose, pH 5.5.
SHIPPING	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Store at≤-70°C, stable for 6 months after receipt. Store at≤-70°C, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
QUALITY	Mol Mass:80.6kDa AP Mol Mass:90-120kDa, reducing conditions.
CONTROL	Purity: Greater than 95% as determined by reducing SDS-PAGE.Endotoxin: Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
BACKGROUND	Glutamate carboxypeptidase2, also known as FOLH1, PSMA, belongs to the M28B subfamily and the peptidase M28 family. It is highly expressed in prostate epithelium and can be detected in urinary bladder, kidney, testis, ovary, fallopian tube, breast, adrenal gland, liver, esophagus, stomach, small intestine, colon and brain (at protein level). PSMA is used as a diagnostic and prognostic indicator of prostate cancer, and as a possible marker for various neurological disorders such as schizophrenia, Alzheimer disease and Huntington disease. It has both folate hydrolase and N-acetylated-alpha-linked-acidic dipeptidase (NAALADase) activity and has a preference for tri-alpha-glutamate peptides. PSMA involves in prostate tumor progression and also exhibits a dipeptidyl-peptidase IV type activity. In vitro, PSMA cleaves Gly-Pro-AMC. PSMA is stable at pH greater than 6.5.
KDa MK R 120 90 60 5DS-PAGE 40 30 20	