

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

Recombinant Human TROP-2

Catalog#:P00490 Derived from Human Cells

DESCRIPTION	Recombinant Human Tumor-associated Calcium Signal Transducer 2 is
	produced by our Mammalian expression system and the target gene encoding
	Thr88-Thr274 is expressed with a 6His tag at the C- terminus.
	Accession#: P09758
	Known as : Tumor-associated calcium signal transducer2; Membrane component
	chromosome 1 marker 1; Cell surface glycoprotein Trop-2; TACSTD2; TROP2
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, 2mM EDTA,
	pH 7.4.
SHIPPING	The product is shipped at ambient temperature.
	Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at <-20°C, though stable at room
	temperature for 3 weeks.
	Reconstituted protein solution can be stored at 4-7°C for 2-7 days.
	Aliquots of reconstituted samples are stable at < -20°C for 3 months.
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting.
	It is not recommended to reconstitute to a concentration less than 100µg/ml.
	Dissolve the lyophilized protein in distilled water.
	Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY	Mol Mass:22.6kDa AP Mol Mass:28-40kDa, reducing conditions.
	Purity : Greater than 95% as determined by reducing SDS-PAGE.
CONTROL	Endotoxin: Less than 0.1 ng/μg (1 EU/μg) as determined by LAL test.
BACKGROUND	Tumor associated calcium signal transducer 2 (TACSTD2, TROP-2) is a type I cell surface glycoprotein that is highly expressed on human carcinomas. It was originally identified as an antigen present on human gastrointestinal tumors and is the second of two members of this family. Human and mouse TROP-2 share 87% amino acid (aa) similarity. TROP-2 is capable of transducing an intracellular
	originally identified as an antigen present on human gastrointestinal tumors and
	is the second of two members of this family. Human and mouse TROP-2 share
	87% amino acid (aa) similarity. TROP-2 is capable of transducing an intracellular
	calcium signal and may play a role in tumor growth. It also has adhesive functions.
	I TUTIONIS.



