

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

Recombinant Human SEMA3A

Catalog#:P01432 Derived from Human Cells

Catalog#.1 01432 Derived Holli Hullian Cens	
DESCRIPTION	Recombinant Human Semaphorin 3A is produced by our Mammalian expression system and the target gene encoding Lys26-Phe546 is expressed with a 6His-Flag tag at the C-terminus. Accession#: Q14563 Known as: Semaphorin III; SEMA3A; SEMAD; Hsema- I; SEMA1; SemD; Semad; coll- 1; Hsema- III; Sema III
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, 5% Trehalose, pH 7.4.
SHIPPING	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at<-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY	Mol Mass:61.3kDa AP Mol Mass:70-75kDa, reducing conditions. Purity: Greater than 95% as determined by reducing SDS-PAGE.
CONTROL	Endotoxin : Less than 0.1 ng/ μ g (1 EU/ μ g) as determined by LAL test.
BACKGROUND	emaphorin-3A is a secreted protein which belongs to the semaphorin family. Semaphorins are a family of secreted and cell-bound signaling molecules defined by the presence of a common 500 aa Sema domain. This protein can function as either a chemorepulsive agent, inhibiting axonal outgrowth, or as a chemoattractive agent, stimulating the growth of apical dendrites. In both cases, the protein is vital for normal neuronal pattern development. Increased expression of this protein is associated with schizophrenia and is seen in a variety of human tumor cell lines.
	KDa MK R 120 90 60 SDS-PAGE