

Tel: 400-968-6088 Fax: 010-56371281

Http://www.solarbio.com

## **Recombinant Mouse IL-21**

Catalog#:P02069 Derived from *E.coli* 

Catalog#.F02009 Delived Holli E.coli	
DESCRIPTION	Recombinant Mouse Interleukin-21 is produced by our E.coli expression system and the target gene encoding Pro25-Ser146 is expressed.  Accession#: Q9ES17  Known as: Interleukin-21; Il21
FORMULATION	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4.
SHIPPING	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature listed below.
STORAGE	Lyophilized protein should be stored at<-20°C, though stable at room temperature for 3 weeks.  Reconstituted protein solution can be stored at 4-7°C for 2-7 days.  Aliquots of reconstituted samples are stable at < -20°C for 3 months.
RECONSTITUTION	Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in distilled water. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.
QUALITY CONTROL	Mol Mass: 14.4kDa AP Mol Mass: 16kDa, reducing conditions.  Purity: Greater than 95% as determined by reducing SDS-PAGE.  Endotoxin: Less than 0.001 ng/ug (0.01 EU/ug) as determined by LAL test.
BACKGROUND	Interleukin-21 also known as IL-21 is a protein that in mouse is encoded by the IL21 gene, belongs to the IL-15/IL-21 family. Interleukin-21 is a cytokine that has potent regulatory effects on cells of the immune system, including natural killer (NK) cells and cytotoxic T cells that can destroy virally infected or cancerous cells. This cytokine induces cell division/proliferation in its target cells.
kDa MK R 120 90	

