

Product Specification

Cat No:	BYA1031
English Name:	17 Amino Acids
Size:	1mL
Storage:	Store at 2-8°C, 2 years.

Component Information:

CAS No.	Name	Standard value	Unit
56-84-8	L-Aspartic Acid	1	mmol/L
72-19-5	L-Threonine	1	mmol/L
56-45-1	Serine	1	mmol/L
56-86-0	L-Glutamic acid	1	mmol/L
147-85-3	L-Proline	1	mmol/L
56-40-6	L-Glycine	1	mmol/L
338-69-2	D-Alanine	1	mmol/L
56-89-3	L-Cystine	1	mmol/L
72-18-4	L-Valine	1	mmol/L
59-51-8	DL-Methionine	1	mmol/L
73-32-5	L-Isoleucine	1	mmol/L
61-90-5	L-Leucin	1	mmol/L

60-18-4	L-Tyrosine	1	mmol/L
63-91-2	L-Phenylalanine	1	mmol/L
56-87-1	L-Lysine	1	mmol/L
71-00-1	L-Histidine	1	mmol/L
74-79-3	L-Arginine	1	mmol/L

References:

Welcome to follow our company scholarship program

Related products:**BYA9010**

Amino acid (20AA)

SDK1010

Amino acid derivatization kit (OPA method, without standards)

SDK1020

Amino acid derivatization kit (PITC method, without standards)

SG8540

L-Glutamic acid

SL8330

L-Lysine

SA9160

L-Alanine