

Anti-TCEAL3 Polyclonal Antibody

Cat: K107688P


Summary:

【Product name】 : Anti-TCEAL3 antibody	【Source】 : Rabbit
【Isotype】 : IgG	【Species reactivity】 : Human Mouse
【Swiss Prot】 : Q969E4/Q5H9L2/Q6IPX3	【Gene ID】 : 85012
【Calculated】 : MW:22kDa	【Observed】 : MW:32kDa
【Purification】 : Affinity purification	
【Tested applications】 : WB IHC	
【Recommended dilution】 : WB 1:5000-1:8000. IHC 1:25-100.	
【WB Positive sample】 : K-562,A431,SH-SY5Y	
【IHC Positive sample】 : Human colorectal cancer	
【Subcellular location】 : Cytoplasm Secreted	
【Immunogen】 : A synthetic peptide of human TCEAL3	
【Storage】 : Shipped at 4°C. Upon delivery aliquot and store at -20°C	


Background:

May be involved in transcriptional regulation.

Verified picture



Western blot analysis with TCEAL3 antibody diluted at 1:6000; Lane: K-562, A431, SH-SY5Y



Immunohistochemistry of paraffin-embedded Human colorectal cancer with TCEAL3 antibody diluted at 1:40