

## Anti-RAC3 Polyclonal Antibody

Cat: K108127P

### Summary:

---

**【Product name】** : Anti-RAC3 antibody

**【Source】** : Rabbit

**【Isotype】** : IgG

**【Species reactivity】** : Human Mouse Rat

**【Swiss Prot】** : P60763

**【Gene ID】** : 5881

**【Calculated】** : MW:21kDa

**【Purification】** : Affinity purification

**【Tested applications】** : IHC

**【Recommended dilution】** : IHC 1:50-200.

**【IHC Positive sample】** : Human breast cancer

**【Subcellular location】** : Cytoplasm

**【Immunogen】** : Recombinant protein of human RAC3


**【Storage】** : Shipped at 4°C. Upon delivery aliquot and store at -20°C

### Background:

---

The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Alternative splicing results in multiple transcript variants.

## Verified picture



Immunohistochemistry of paraffin-embedded  
Human breast cancer with RAC3 antibody  
diluted at 1:100