

Hexadecylpyridinium chloride monohydrate

Cat: C9890

Specification: 25g

Storage: Store at room temperature, and it is valid for 2 years.

Product Information

CAS: 6004-24-6

English name: Hexadecylpyridinium chloride monohydrate

Molecular Formula: C₂₁H₃₈ClN·H₂O

Molecular Weight: 358

Appearance (Character): White crystalline powder

Purity: ≥ 99%

Solubility: 179 mg/mL in Water

Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.