

## Didodecyldimethylammonium Bromide

**Cat:** D7370

**Storage:** Store at room temperature, avoid light, and it is valid for 2 years.

### Product Information

**CAS:** 3282-73-3

**English name:** Didodecyldimethylammonium bromide

**Molecular Formula:** C<sub>26</sub>H<sub>56</sub>BrN

**Molecular Weight:** 462.63

**Appearance (Character):** White crystalline powder

**Purity:** ≥ 98%

**Solubility:** dissolved in tritiated chloroform

### Introduction

Phase transfer catalysts are widely used as cracking reagents in pharmaceuticals and clinical diagnostic kits; It is also used as an emulsifier, bactericidal agent, disinfectant, and antistatic agent. such as bioengineering and electronics.

### Note

1. When taking it, it must be dry, avoiding frequent dissolution and freeze-drying.
2. This product is only for scientific research and cannot be used for diagnosis and treatment.