

## Glycogen from oyster TypeII

**Cat:** G8050

**Specification:** 1g

**Storage:** Store at 2-8°C, and it is valid for 2 years.

### Product Information

**CAS:** 9005-79-2

**English name:** Glycogen from oyster TypeII

**Appearance (Character):** White powder

**Source:** From Oyster

**Type:** Substrate of Type II enzyme

**Solubility:** 50mg/ml in H<sub>2</sub>O

### Introduction

Glycogen is an intermediate product of glycolysis and a high-energy phosphate, serving as a form of energy storage for humans, animals, fungi, and bacteria. It can be used for sugar phosphorylation.

### Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
3. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
4. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.