

Pituitary Luteinizing Hormone

Cat: L8040

Specification: 200 units

Storage: Keep in a cool, dark place, away from light and sealed, and it is valid for 2 years.

Product Information

English name: Pituitary Luteinizing Hormone

Appearance (Character): White or off-white freeze-dried blocks or powder.

Pharmacological Action: In synergy with pituitary follicle-stimulating hormone, this product can promote the final maturation of ovarian follicles and induce the formation of mature follicles and corpora lutea.

Indications: Used for the treatment of ovulation delay, ovarian cysts, habitual abortion, etc.

Dosage and Administration: Intramuscular Injection: Single dose for horses is 200-300 units; for cattle, 100-200 units. Before use, dilute with 2-5ml of sterile normal saline.

Adverse Reactions: May cause allergic reactions.

Note: When treating ovarian cysts, the dosage should be doubled.