

Olive oil

Cat: O8320

Specification: 500mL

Storage: Store at room temperature, and it is valid for 2 years.

Product Information

CAS: 8001-25-0

English name: Olive oil

Appearance (Character): Yellowish green or yellowish brown oily liquid

Acid value: ≤ 4.0 (mgKOH/g)

Density: 0.910-0.915(20°C)/(g/ml)

Iodine value: 79.0-88.0

Refractive index: 1.468-1.470

Grade: CP

Introduction

It is a biochemical reagent that can be used as a biomaterial or organic compound for life science-related research.

Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
3. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
4. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.