

PMSG

Cat: P9970

Specification:5*1000IU

Storage: Store at 4°C, avoid light, and it is valid for 2 years.

Product Information English name: PMSG

Appearance (Character): White solid

Introduction

Function: Hormone drugs. It has the effect of promoting follicle stimulating hormone and luteinizing hormone. It can promote the development and maturation of follicles in female livestock, causing them to enter estrus; Promote the ovulation of mature follicles and even superovulation. Increase the secretion of androgens in male animals and enhance sexual excitement.

Purpose: 1 Superovulation during embryo transfer

2. Stimulating estrus and promoting follicular development in female animals

Usage: Before use, dilute it with 2-5 mL of sterile saline for subcutaneous or intramuscular injection.

The dosage is for reference only:

Aphrodisiac: Horse and cow: 1000-2000 units;

Sheep: 100-500 units; Pig: 200-800 units; Canine: 25-200 units; Cat: 25-100 units;

Rabbit and mink: 30-50 units;

Superovulation: cows: 2000-4000 units;

Female sheep: 600-1000 units;

Note

- 1. Long-term use should be avoided to prevent the production of antibodies and inhibition of pituitary gonadotropic function.
- 2. The solution of this product is extremely unstable and heat-sensitive, so it should be used up within a short period of time.
- 3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
- 4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
- 5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.