

Polyethylene-polypropylene glycol

Cat: S7070

Specification: 100g

Storage: Store at room temperature, and it is valid for 2 years.

Product Information

CAS: 9003-11-6

English name: Polyethylene-polypropylene glycol

Appearance (Character): White powder or granules

Solubility: 100mg/ml in Water

Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.