

## Surfactant S7

**Cat:** S7420

**Specification:** 10mL /50mL /100mL

**Storage:** Store at room temperature, and it is valid for 2 years.

### Product Information

**English name:** Surfactant S7

**Appearance (Character):** Colorless liquid

**Category:** Anionic

**pH:** 6.0~8.0

**Viscosity:** 26-36 mPas

### Introduction

Surfactant S7 is a potent anionic surfactant that can be used as a surfactant in rapid diagnostic test strips. It is specifically applied as a surfactant for colloidal gold.

### Applicable Samples:

Suitable for a full range of products, with particularly better performance in viscous oral samples.

### Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.