

## Surfactant S17

**Cat:** S7440

**Specification:** 25g /500g

**Storage:** Store at room temperature, and it is valid for 2 years.

### Product Information

**English name:** Surfactant S17

**Appearance (Character):** White or light yellow

**Category:** Non-ionic

**Molecular Weight:** 1148

**pH Value:** 6-8 at 10wt/wt%

**Density:** 1.03~1.04 mg/ml at 25°C (77°F)

**Moisture Content:** ≤0.50%

**HLB Value:** 14

### Introduction

Surfactant S17, a non-ionic surfactant, is an excellent emulsifier, solubilizer, and dispersant.

### Applicable Samples:

Suitable for a full range of products, with particularly better performance in urine samples.

### Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.