

## Tryptone Soy Broth

**Cat:** LA0110

**Specification:** 250g

**Storage:** Store at room temperature, and it is valid for 3 years.

### Product Information

**English name:** Tryptone Soy Broth

**Alias:** TSB; Tryptone Soy Broth TSB Medium

**Appearance (Character):** Off-White powder

**Grade:** BR

### Introduction

For MPN determination of Staphylococcus aureus, enrichment culture.

### Composition: g/L

|                                 |     |
|---------------------------------|-----|
| Tryptone:                       | 17  |
| Phytopeptone:                   | 3   |
| Sodium chloride:                | 5   |
| Dipotassium hydrogen phosphate: | 2.5 |
| Glucose:                        | 2.5 |

### Instructions for use(*for reference only*):

Weigh 30g of this product, add it to 1000ml distilled water, dissolve it by boiling, divide into aliquots, autoclaving at 121 ° C for 15 minutes, and set to use(The concentration of sodium chloride can be adjusted as needed; Do not use this product if clumps occur; After sterilization, the product should be operated in a sterile environment).

### Note

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
3. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
4. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.