

Hoagland Nutrient Solution (Normal, Calcium Nitrate Free)

Cat: LA2061

Specification: 250g

Storage: Store at 2-8°C, and it is valid for 2 years.

Product Information:

Chinese name: 霍格兰营养液

Composition (mg/L) :

Potassium sulfate	607
Ammonium dihydrogen phosphate	115
Magnesium sulfate	493
Ferric sodium salt	20
Ferrous sulfate	15
Boric acid	2.86
Borax sand	4.5
Manganese sulfate	2.13
Copper sulfate	0.05
Zinc sulfate	0.22
Ammonium sulfate	0.02

Introduction:

Used as a nutrient solution for soilless cultivation of plants.

Instructions for use: (for reference only)

Weigh 1.26g of this product, and take 0.945g of calcium nitrate, dissolve in 1000ml distilled water by heating, repack, autoclast sterilization at 115°C for 20 minutes, and set to use.

Note:

1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
2. Once dissolved, please store the solution in separate containers to avoid product degradation caused by repeated freezing and thawing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.