

## **BMGY/BMMY Medium (ready-to-use)**

Cat No.: LA5252

Package: 500mL

**Storage:** Store at  $4-8^{\circ}$ C for one year.

# **English name:** BMGY/BMMY Medium **Product brief introduction**

It is used to induce the expression of Pichia Pastoris, containing Biotin and non-amino yeast nitrogen (YNB).

### Product formula and usage:

Product ingredient	Concentration
Yeast Extract Powder	1%
Peptone	2%
Amino free yeast nitrogen Source (YNB)	1.34%
Biotin (D-biotin)	0.4mg/L
Phosphate buffer system	100mM

It can be used to induce the expression of Pichia pastoris by adding 1% glycerol or 0.5% methanol. Usage: Disinfect the bottle cap with 75% alcohol, and use it directly after opening the plastic cap under sterile conditions.

#### This reagent is only used in scientific research.

#### **Related products**

LA0223 YPD Liquid Medium (Ready-to-use)
Y8063 YPDA Liquid Medium (Ready-to-use)
LA0661 MD Liquid Medium Minimal Dextrose Medium
LA5252 BMGY/BMMY Medium (ready-to-use)

Note: For more information about this product, please refer to Solarbio website.

