

Pseudomonas CFC Selective Agar additive

Cat: LA6061

Specification: 1ml×20 pieces/box

Storage: Store at 2-8°C, and it is valid for 1 year.

Product Information:

Main ingredients:

cetane trimethyl ammonium bromide, 1.0mg; Sodium fusidate, 1.0mg; Cephalosporin II, 5.0mg;
Distilled water, 1.0ml

Performance indicators:

1. Physical and chemical properties: colorless or slightly yellow, with possible slight turbidity but not affecting use.
2. Bacteriostasis test: Qualified.

Instructions for use: (*for reference only*)

In a sterile environment, use a grinding wheel to cut the bottleneck directly below the blue dot on the ampoule bottle, and open the ampoule bottle with your back facing away from the blue dot while wearing gloves. When the culture medium has cooled to 50-55°C, aseptically add an appropriate amount of this product to the culture medium using a pipette, gently mix it, and set it aside.

Note:

1. This product is a supporting additive for culture medium, please use it together.
2. This product is only used for laboratory microbiological testing.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.