

Thioglycollate Broth Medium

Cat: LA8740

Specification: 250g

Storage: Store at room temperature, and it is valid for 3 years.

Product Information:

English name: Thioglycollate Broth Medium

Appearance (Character): Light yellow powder

pH: 6.9-7.3

Grade: BR

Composition: g/L

Casein peptone:	15.0
Yeast infusion powder:	5.0
L-cystine:	0.5
Sodium thioglycolate:	0.5
Sodium chloride:	2.5
Glucose:	5.5
pH	7.1±0.2, 25°C

Instructions for use: (*for reference only*)

Weigh 29g of this product into 1L of distilled water or deionized water, boil and dissolve, divide into packages, sterilize at 121°C for 15 minutes or at 115°C for 30 minutes, and cool quickly after sterilization.

Note:

1. This product is not suitable for intraperitoneal injection in mice to induce macrophage formation.
2. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
3. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
4. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
5. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.