

## Potassium Tellurite Egg-Yolk Reagent

**Cat:** RL100109

**Specification:** 6\*25ml

**Storage:** Store at 2-8°C, valid for 2 years

### Introduction:

#### I. Product use:

Covered liquid medium for anaerobic culture.

#### II. Composition: g/L

Formula (per piece) Content

30% vitelline saline: 20.83mL

1% potassium tellurite solution: 4.17mL

#### III. Instructions for use: (for reference only)

Each one was added in 475mLBaird-Parker Agar base medium.

#### IV. Note:

Operate in a clean environment to avoid contamination of the culture medium.

#### V. Waste disposal:

After testing, the contaminated items are placed under high-pressure sterilization at 121°C for 30 minutes.