

Soybean protein isolate

Cat: S9510

Specification: 250g

Storage: Store at room temperature, and it is valid for 3 years.

Product Information

Appearance (Character): Milky white or creamy yellow powder

Protein: ≥85% Moisture: ≤7%

Ash content: ≤15% pH value: 5.5-7.5

Grade: BR

Introduction

Soy Protein Isolate contains over 90% protein, with nearly 20 types of amino acids, including various essential amino acids. It is nutritious and suitable for biological culture research.

Note: This product is not fully soluble and can be used for fermentation.

Note

- 1. Unless otherwise specified, the biochemical reagents produced by our company are generally non-sterile packaged. If they are to be used for cell experiments, please conduct pretreatment in advance.
- 2. The product information is for reference only. If you have any questions, please call 400-968-6088 for consultation.
- 3. The products are all for scientific research use only. Do not use it for medical, clinical diagnosis or treatment, food and cosmetics, etc. Do not store them in ordinary residential areas.
- 4. For your safety and health, please wear laboratory clothes, disposable gloves and masks to operate.